

COLLECTION TUBE CODE MASTER LISTING

BLOOD Collection Tube Code	Name	NON BLOOD Specimen Collection Code	Name	Tolerance Testing Collection Tube Code	Name
1L9	Ref Lavender WB, Aux1	1OB	Fecal Occult Blood Card, 1st Specimen	T1C	Tolerance Sample 1
2L9	Ref Lavender WB, Aux2				
A1S	SST,Immulate1				
A2S	SST,Immulate 2				
A3S	SST,Immulate 3				
AAD	Yellow ACD, Solution A				
ABD	Yellow ACD Solution B				
ACD	Yellow ACD, Solution A or B				
ALS	SST - Immulate Testing				
B1U	Blue, Na Citrate, Second Tube				
B2U	Blue, Na Citrate, Third Tube				
BLU	Blue, Na Citrate				
BNS	Blue, Na Citrate, No Spin				
COL	Lavender,EDTA				
CTC	CellSearch CellSave Tube, Room Temp	2OB	Fecal Occult Blood Card, 2nd Specimen	T2C	Tolerance Sample 2
DBP	Dark Blue, w/EDTA	3OB	Fecal Occult Blood Card, 3rd Specimen	T3C	Tolerance Sample 3
DER	Dark Blue w/EDTA, Refrigerated	APC	AP container	T4C	Tolerance Sample 4
DRB	Dark Blue, No Additive	BFL	Body Fluid, Lavender Tube	T5C	Tolerance Sample 5
GRN	Green, Lithium Heparin,WB	CF0	CSF	T6C	Tolerance Sample 6
GYT	Gray	ELL	Specimen to Elliot	T7C	Tolerance Sample 7
ICA	SST for Ionized Calcium	FLC	Body Fluid Collection Container	T8C	Tolerance Sample 8
L1V	Lavender, EDTA	GAS	Resp Blood Gas	U10	Urine Tolerance Sample 1
L2V	Lavender, EDTA	HBT	Breath Tek-UBit-UBT Kit	U20	Urine Tolerance Sample 2
LAV	Lavender, EDTA	ONC	Non Sterile Specimen Container	U30	Urine Tolerance Sample 3
LBB	Blood Bank Lavender	OPF	O&P, 10%Formalin, &PVA	U40	Urine Tolerance Sample 4
		RGL	Reagan-Lowe Transport Tube-Bordetella	U50	Urine Tolerance Sample 5
		SPC	Specimen Container	U60	Urine Tolerance Sample 6
		SUR	SurePath	U70	Urine Tolerance Sample 7
		SWB	Swab		
		SYR	Bld Gas Syringe		
		THN	Thin Prep Vial		
		UA0	Urine Specimen		
		VCM	Viral-Chlam-Mycoplasma transport media		

COLLECTION TUBE CODE MASTER LISTING, CONT'D

BLOOD Collection Tube Code	Name	NON BLOOD Specimen Collection Code	Name	Tolerance Testing Collection Tube Code	Name
LCD	Lavender, EDTA, Cord Blood	VTM	Viral Transport Media		
LEG	Legal Specimen	XCF	CSF for Micro, Last Tube		
LRR	Lavender, EDTA, Refrigerated	XGE	Swab (Gel)		
MIB	Microbilirubin, Capillary	XMT	Swab (Mini-Tip)		
MSL	Misc Specimen Container - Call lab for instructions	XNG	Swab (Non-gel)		
N1G	Green, SODIUM Heparin	XRC	Sterile Container, Refrigerated		
N2G	Green, SODIUM Heparin				
NAG	Green, SODIUM Heparin				
NH3	PST, on Ice for Ammonia	XST	Sterile Container		
NSP	Green, NA Hep, No Spin	XXF	Sterile Container, Frozen		
NUL	No Specimen				
P1T	PST, Second Tube Needed				
P2T	PST, Third Tube Needed				
PLV	Lavender, Spun				
PPT	White Top, PPT, EDTA				
PST	Plasma Separator Tube				
R1D	Red, Second Tube Needed				
R2D	Red, Third Tube Needed				
RBB	Blood Bank Red				
RED	Red				
RLB	Blood Bank Red and Lavender				
S1T	SST, Second Tube Needed				
S2T	SST, Third Tube Needed				
SST	SST				
VAP	SST				

COLLECTION TUBE CODE MASTER LISTING, CONT'D

BLOOD Collection Tube Code	Name		NON BLOOD Specimen Collection Code	Name		Tolerance Testing Collection Tube Code	Name
XBC	Blood Culture						
XFU	Blood Culture, Fungal						
XPS	SPS, Sterile Yellow top, Blood Collection Tube						
Z19	Blue, Na Citrate, Second Tube Needed, Ref Lab Test						
Z29	Blue, Na Citrate, Third Tube Needed, Ref Lab Test						
Z39	Ref Blue Na Cit, Fourth tube						
ZAW	Lavender, EDTA, WB, Reference Lab Test, Room Temp						
ZB9	Blue, Na Citrate, Reference Lab Test						
ZL9	Lavender, EDTA, WB, Reference Lab Test						
ZR9	Reference Lab Red Top						
ZWF	Lavender, EDTA, WB, Reference Lab Test, Frozen						

COLLECTION TUBE CODE MASTER LISTING, CONT'D